



Digital Egypt Pioneers Initiative Pneumonia & Respiratory Disease Severity Scoring

Project Report

Abdellatif Elbatrawy

Omar Essam

Ali Mohamed

Amr Ghoniem

Ibrahim Hanafy

Ramez Farouk

Supervisor

Eng. Abdelrahman El Mashtoly

Contents

1	Executive Summary	2
1.1	Overview	2
1.2	Objective	2
1.3	Technical Approach	2
1.4	Key Results & Impact	2
1.5	Deployment	2
2	Methodology	3
2.1	Data Strategy	3
2.2	Preprocessing & Augmentation	3
2.3	Model Evolution	3
2.3.1	Phase A: The Baseline (Custom CNN)	4
2.3.2	Phase B: ResNet50V2	4
2.3.3	Phase C: The Ensemble Method	5
2.3.4	Phase D: DenseNet121	5
2.4	Training Pipeline: Sequential Transfer Learning	6
3	System Design	7
3.1	Technology Stack	7
3.2	Application Workflow	7
4	Results & Performance	8
4.1	Quantitative Metrics	8
4.2	Analysis of Results	9
4.3	Visualizations	9
5	Conclusion & Future Scope	9
5.1	Conclusion	9
5.2	Future Work	10
	Appendix: User Manual	11

1 Executive Summary

1.1 Overview

Respiratory diseases, particularly pneumonia, remain a critical public health challenge globally and in Egypt. Pneumonia is a leading cause of mortality among children and the elderly, claiming over 2.5 million lives annually worldwide. The diagnostic process relies heavily on Chest X-rays (CXR), which require expert interpretation by radiologists. However, the shortage of radiologists, coupled with high patient volumes in public hospitals, often leads to diagnostic delays and fatigue-induced errors. Studies suggest that missed pneumonia diagnoses (false negatives) can occur in up to 20% of cases, leading to severe complications or death.

1.2 Objective

The primary goal of this project, developed under the Digital Egypt Pioneers Initiative (DEPI), was to engineer an automated, high-sensitivity Artificial Intelligence system capable of acting as a "Universal Radiologist." Unlike traditional AI models that are often biased toward specific demographics (e.g., adults only), our objective was to build a robust system that performs reliably across both pediatric and adult populations.

1.3 Technical Approach

Our solution utilizes DenseNet121 Deep Learning architecture, selected for its superior ability to capture fine-grained medical textures. To achieve true "Universal" applicability, we implemented a **Multi-Stage Sequential Transfer Learning** pipeline, training the model on three of the world's largest public medical datasets:

1. **RSNA Pneumonia Challenge:** Established the baseline feature extraction on pediatric data.
2. **NIH Chest X-ray:** Adapted the model to adult anatomy to prevent demographic bias.
3. **CheXpert:** Finalized the model by fine-tuning on this massive dataset (200,000+ images) to ensure robustness against varied imaging conditions, angles, and hospital equipment types.

1.4 Key Results & Impact

The final "Universal Model" achieves diagnostic excellence with an **AUC of 0.89**. Crucially, we calibrated the system to prioritize patient safety over raw accuracy. By optimizing the decision threshold to **0.15**, we achieved a **Sensitivity (Recall) of ~91%** across both adult and pediatric test sets. This ensures that 9 out of 10 potential pneumonia cases are flagged for doctor review, significantly reducing the risk of missed diagnoses.

1.5 Deployment

The system is deployed as a lightweight, online Python application that runs on standard laptops. This design allows it to be used in remote or resource-limited clinics in Egypt, providing immediate diagnostic support without requiring expensive cloud infrastructure. This project

represents a scalable, digital transformation solution that directly addresses a vital healthcare need.

2 Methodology

2.1 Data Strategy

We utilized a multi-source data approach to ensure robustness:

- **RSNA Pneumonia Challenge:** 26,000 images (Pediatric focus). Used for initial feature learning.
- **NIH Chest X-ray:** 112,000 images (Adult focus). Used for domain adaptation.
- **CheXpert (Stanford):** 224,000 images. Used for robustness and finalization.

2.2 Preprocessing & Augmentation

All images undergo a standardized pipeline before inference to ensure consistency.

- **Resizing:** Upgraded from standard 224×224 to 320×320 to capture fine-grained lung textures.
- **Normalization:** Pixel intensity scaled to $[0, 1]$.
- **Augmentation:** To prevent overfitting, we applied random transformations during training.

```
1 train_datagen = ImageDataGenerator(  
2     rescale=1./255,  
3     rotation_range=15,  
4     zoom_range=0.2,  
5     horizontal_flip=True,  
6     fill_mode='nearest'  
7 )
```

Listing 1: Data Augmentation Configuration




Figure 1: Augmentation on the Pneumonia case

2.3 Model Evolution

We did not simply choose a model; we engineered the optimal solution through four distinct experimental phases.

2.3.1 Phase A: The Baseline (Custom CNN)

We began by training a custom Convolutional Neural Network from scratch.

- **Result:** The model plateaued at 76% Accuracy.
- **Analysis:** The network lacked the depth to extract complex features like lung consolidation.

```
1 model = Sequential([
2     Conv2D(32, (3,3), activation='relu', input_shape=(224, 224, 3)),
3     MaxPooling2D(2,2),
4     Conv2D(64, (3,3), activation='relu'),
5     MaxPooling2D(2,2),
6     Flatten(),
7     Dense(128, activation='relu'),
8     Dropout(0.5),
9     Dense(1, activation='sigmoid')
10 ])
```

Listing 2: Custom CNN Architecture

2.3.2 Phase B: ResNet50V2

Understanding the need for deeper feature extraction, we implemented Transfer Learning using ResNet50V2, a state-of-the-art residual network pre-trained on ImageNet.

- **Optimization:** We utilized a Rescaling layer to normalize inputs to the specific range $[-1, 1]$ required by ResNet.
- **Result:** Performance improved to 78.5% Accuracy. While better than the baseline, ResNet50 struggled to distinguish subtle pneumonia cases from other lung noise.

```
1 inputs = tf.keras.Input(shape=(320, 320, 3))
2 x = Rescaling(scale=2.0, offset=-1.0)(inputs) # Preprocessing [-1, 1]
3
4 # Load Backbone (Frozen)
5 base_model = ResNet50V2(weights='imagenet', include_top=False,
6     input_tensor=x)
7 base_model.trainable = False
8
9 # Custom Head
10 x = base_model.output
11 x = GlobalAveragePooling2D()(x)
12 x = Dense(256, activation='relu')(x)
13 x = Dropout(0.5)(x)
14 predictions = Dense(1, activation='sigmoid')(x)
15 model = Model(inputs=inputs, outputs=predictions)
```

Listing 3: ResNet50V2 Implementation

2.3.3 Phase C: The Ensemble Method

We attempted to boost performance by creating a Weighted Ensemble of the DenseNet121 (60% weight) and ResNet50V2 (40% weight).

- **Hypothesis:** Combining two architectures would cancel out individual errors.
- **Result:** The ensemble achieved 85.02% Accuracy, slightly lower than the single DenseNet. We discarded this approach to prioritize inference speed.

```
1 # 1. Generate Predictions from both models
2 p1 = model_dense.predict(test_generator)
3 p2 = model_resnet.predict(test_generator)
4
5 # 2. Apply Weighted Average (60% DenseNet, 40% ResNet)
6 final_ensemble_preds = (p1 * 0.6) + (p2 * 0.4)
```

Listing 4: Ensemble Logic

2.3.4 Phase D: DenseNet121

We selected DenseNet121 as our final architecture due to its superior "feature reuse" mechanism. By connecting every layer to every subsequent layer, DenseNet preserves low-level textural information (edges, haziness) critical for medical imaging.

- **Key Innovation (Stability Fix):** We replaced the standard Sigmoid activation with a Linear (Logits) output to prevent numerical instability (NaN errors) during training.
- **Resolution Upgrade:** We utilized Progressive Resizing, scaling the input from 224×224 to 320×320 to capture fine-grained details.
- **Result:** This architecture achieved 92% Accuracy on Adults and $\sim 90\%$ Sensitivity, significantly outperforming all previous attempts.

```
1 base_model = DenseNet121(weights='imagenet', include_top=False,
2   input_shape=(320, 320, 3))
3 base_model.trainable = False
4
5 x = base_model.output
6 x = GlobalAveragePooling2D()(x)
7 x = BatchNormalization()(x) # Stabilize inputs
8 x = Dense(512, activation='relu')(x)
9 x = Dropout(0.5)(x) # Prevent overfitting
10 x = Dense(256, activation='relu')(x)
11 x = Dropout(0.4)(x)
12
13 predictions = Dense(1, activation=None)(x) # Logits for Stability
14 model = Model(inputs=base_model.input, outputs=predictions)
```

Listing 5: Final DenseNet121 Architecture

2.4 Training Pipeline: Sequential Transfer Learning

Standard transfer learning typically involves training a model on one dataset and stopping. However, given the demographic disparities in pneumonia presentation (pediatric vs. adult), a single-stage training approach proved insufficient. To solve this, we engineered a Multi-Stage Sequential Transfer Learning pipeline.

Stage 1: Feature Extraction

Objective: Adapt ImageNet weights to medical imaging without "shocking" pre-trained extractors.

Dataset: RSNA Pneumonia Detection Challenge (Pediatric focus).

Technique: Froze backbone, trained custom head for 5 epochs.

Outcome: Baseline accuracy of $\sim 76\%$.

Stage 2: Deep Fine-Tuning

Objective: Learn subtle disease textures (opacity, consolidation) specific to children.

Technique: Unfroze all layers. Used SGD with reduced learning rate ($1e^{-4}$) and Gradient Clipping.

Outcome: Accuracy improved to $\sim 85\%$ on pediatric test set.

Stage 3: Domain Adaptation

Objective: Adapt the model to adult anatomy while retaining pediatric knowledge.

Dataset: NIH Chest X-ray Dataset (Adult focus).

Technique: Fine-tuned on adult data using micro-learning rate ($1e^{-5}$).

Outcome: "Universal" model achieved 92.4% Accuracy on Adults while maintaining $\sim 80\%$ Accuracy on Children.

Stage 4: Robustness Finalization

Objective: Stress-test on a massive, noisy, real-world dataset.

Dataset: CheXpert (Stanford Hospital).

Outcome: While this stage improved general robustness, our evaluation showed that the Stage 3 Universal Model offered the best balance of sensitivity.




Figure 2: Sequential Transfer Learning Pipeline Diagram

3 System Design

3.1 Technology Stack

The system was built using a modular Python ecosystem designed for portability and efficiency on standard hardware.

- **Deep Learning Core:** TensorFlow 2.16 / Keras 3.0.
- **Data Processing:** NumPy, Pandas, Pillow (PIL).
- **Visualization:** Matplotlib, Seaborn.
- **Deployment Environment:** Local Python Runtime (Web Application).
- **Hardware Used:** Kaggle NVIDIA Tesla T4 x2 (Distributed Strategy).

3.2 Application Workflow

We developed a standalone website that encapsulates the entire diagnostic logic. The workflow ensures that every input image undergoes the exact same transformation as the training data.

The Logic Flow:

1. **Input:** User provides a local file path to a Chest X-ray (JPEG/PNG).
2. **Preprocessing Engine:**
 - **Conversion:** Grayscale images are converted to RGB (3-channel).
 - **Resizing:** Image is scaled to 320×320 pixels.
 - **Normalization:** Pixel intensity is rescaled to the $[0, 1]$ range.

3. **Model Inference:** The image is passed through the DenseNet121 architecture.
4. **Probabilistic Output:** The model outputs a raw "Logit" score, which is converted via a Sigmoid Activation function into a probability percentage.
5. **Safety Logic (Thresholding):** The probability is compared against our calibrated safety threshold (0.15).
6. **Result:** The system returns a binary diagnosis (Normal or Pneumonia) along with a confidence score.




Figure 3: Enter Caption

4 Results & Performance

4.1 Quantitative Metrics

To validate the "Universal" capability and real-world robustness of our solution, we evaluated the final model on three distinct test sets:

1. **Pediatric (Child):** Unseen data from the Kaggle/RSNA dataset.
2. **Adult (General):** The held-out validation set from the NIH dataset.
3. **External Stress Test (CheXpert):** A completely unseen dataset from Stanford Hospital.

Table 1: Final Model Performance by Demographic


Metric	Pediatric	Adult (NIH)	CheXpert	Meaning
Dataset Source	Kaggle Test	NIH Val	Stanford	Diverse Testing
Sample Size	5,283	5,337	10,000	Significant Validation
AUC Score	0.897	0.861	0.852	High Intelligence
Sensitivity	84.36%	91.67%	78.51%	High Safety
Specificity	78.21%	65.57%	77.16%	Conservative

4.2 Analysis of Results

As shown in Table 2, the Universal DenseNet121 model demonstrates exceptional robustness across age groups.

- **Universal Intelligence (AUC):** The Area Under the Curve (AUC) remains consistently high (> 0.85) across all three datasets. This proves that the model has learned the fundamental visual features of pneumonia independent of patient age or hospital equipment.
- **Safety Calibration:** While the model is capable of 92.4% accuracy on adults (at a standard 0.50 threshold), we deployed it with a safety threshold of **0.15**. This trade-off lowered the specific accuracy but boosted **Sensitivity to $\sim 91\%$** .
- **The CheXpert Test:** The model achieved 0.85 AUC on CheXpert without ever training on it. This confirms the model is not overfitting and is robust enough for deployment.
- **Conclusion:** In a clinical triage setting, missing a pneumonia case (False Negative) is dangerous, while double-checking a healthy patient (False Positive) is merely cautious. Therefore, the high sensitivity values confirm the system’s suitability as a reliable screening tool.

4.3 Visualizations



5 Conclusion & Future Scope

5.1 Conclusion

We successfully engineered a **Universal AI Radiologist** that bridges the gap between pediatric and adult diagnostics. By moving beyond simple model training and implementing a sophisticated **Multi-Stage Sequential Transfer Learning** pipeline, we overcame the limitations of demographic bias. The final system acts as a highly sensitive “second pair of eyes,” potentially saving lives by prioritizing high-risk X-rays.

5.2 Future Work

- **Explainability:** Integrate Grad-CAM heatmaps to visually highlight infection areas.
- **Web Deployment:** Migrate the CLI tool to a FastAPI or Streamlit web dashboard.
- **Severity Grading:** Map confidence scores to clinical severity levels (Mild, Moderate, Severe).

Appendix: User Manual

6.1 System Overview

The **Universal AI Radiologist** is deployed as a fully hosted web application. This eliminates the need for local installation, allowing medical professionals to access the diagnostic tool instantly from any device with a web browser.

- **Platform:** Accessible via a public GitHub/Web URL.
- **Interface:** Simple, drag-and-drop web interface.
- **Compatibility:** Works on Desktop, Tablet, and Mobile browsers.

6.2 How to Use

Step 1: Access the Application

Open your web browser and navigate to the project link:

<https://pneumonia-classifier-xum2.onrender.com/>

Step 2: Upload X-ray

- Click the **"Browse Files"** button or drag an X-ray image directly onto the upload area.
- Supported formats: .jpg, .jpeg, .png.

Step 3: View Diagnosis

Once the image is uploaded, the AI model processes it in the cloud. The result will appear instantly on the screen, showing:

- **Diagnosis:** **PNEUMONIA** or **NORMAL**.
- **Confidence Score:** A percentage indicating the model's certainty (e.g., 94.5%).